


# FarrowWrap™ AD OTS and Custom Order Form

Business: \_\_\_\_\_ Patient: \_\_\_\_\_ PO: \_\_\_\_\_  
 Person Measuring: \_\_\_\_\_ Phone/Fax: \_\_\_\_\_ / \_\_\_\_\_ Date: \_\_\_\_\_

**?** This field indicates a measurement necessary to fit a patient for an OTS (off-the-shelf) component. If the patient will not fit into an OTS component, then the additional measurements are required for a custom component. Mark the number of components desired, OTS or custom, in the appropriate field. You may mix and match custom with OTS (ie OTS footpiece & custom legpiece.) If you know what size OTS garment you want, it is not necessary to fill in the measurements.



Please note that “\_” in the SKU below (ie FW\_-O-LR) stands for either “CL” for Classic or “LT” for LITE. The skus may then be cross-referenced with our price sheet to obtain the garment price.

OTS Components		Legpiece									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(C) Mid-Calf		36–43 cm		42–50 cm		48–58 cm		53–63 cm		58–68 cm	
(B) Ankle		21–25 cm		25–30 cm		30–36 cm		36–42 cm		42–50 cm	
(A-D) Regular	FW_-O-LR	35 cm		37 cm		39 cm		41 cm		41 cm	
(B-D)		28 cm		30 cm		31 cm		33 cm		33 cm	
(A-D) Tall	FW_-O-LT	40 cm		42 cm		44 cm		46 cm		46 cm	
(B-D)		33 cm		35 cm		36 cm		38 cm		38 cm	
		Footpiece									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(A <sup>1</sup> ) Mid-Foot		22–24 cm		25–27 cm		28–30 cm		31–34 cm		35–40 cm	
(X) Regular	FW_-O-FR	16.5 cm		18 cm		19.5 cm		21 cm		23 cm	
(X) Long	FW_-O-FL	18.5 cm		20 cm		21.5 cm		23 cm		25 cm	

Options	
Item	√
Classic (Tan)	
LITE (Tan)	
Use Nonstandard Liner Option (see Liner Form)	
When a footpiece & legpiece are purchased together, 1 pair of <b>standard</b> liners are included automatically free. Other liners may cost extra. See directions, liner form, & pricing for details.	

Custom Components	Legpiece	FW_-C-L		Footpiece	FW_-C-F		Extend footpiece length to base of toes
Quantity	Right:		Left:	Right:		Left:	L   R   BL:

Please note that you may substitute a FarrowHybrid™ AD Stocking for a footpiece. See FarrowHybrid™ order form for additional details.

Other: \_\_\_\_\_